

# vPad-IN

## Infant Incubator & Radiant Warmer Testing System

Datrend Systems next generation Incubator Tester, vPad-IN, is based on our revolutionary Vision Pad technology. Utilizing the Vision Pad platform, parameters are sensed inside the incubator environment and results are transmitted wirelessly to a 10" Android tablet for analysis and display. Using wireless control, test parameters and test protocols can be modified without disturbing the equilibrium of the heated environment, thereby minimizing test times.

### Key Features:

- Compatible with closed, forced air convection incubators, and open, radiant infant warmers (optional)
- Simultaneous measurement of humidity, airflow, sound, and 6 different temperatures
- The only Tester available to test in compliance with IEC 60101-2-19 and IEC 60601-2-21
- All sensors are independently movable, to be placed in positions required by the Standard
- Sound level measurements may be taken outside the incubator for external alarm test
- Results can be merged with vPad Safety Test results and imported into commonly-used Equipment Management Systems
- Optional precision controlled heater permits testing of body temperature sensor accuracy

*vPad-IN was designed  
to help protect our  
most precious resource  
in life — the beginning  
of a new life itself!*

Innovation by design



## Incubtest – Performance Specifications

### Temperature:

Measurement Range: 0° C to 50° C

Accuracy: T1+/- 0.05° C

Resolution: 0.01° C

Convection: 5 “puck” sensors (T1-T5) - radiant warmer

Air Convection: 5 sensors (alternate “pos” T1 - T5) - incubator  
1 Mattress Temperature Sensor

Skin Sensor<sup>2</sup>: Sensor “oven” with thermal pad contact

28° C to 40° C in 2° C increments

Resolution: 0.01° C

Accuracy: +/- 0.05° C

### Airflow:

Measurement Range: 0.1 to 1 m/s

Accuracy: +/- 0.1 m/s at 50% RH +/- 15%

### Sound Level:

Measurement Range: 30 dbA to 90 dbA, Type III

Accuracy: +/- 3 dbA

Resolution: 0.1 dbA

### Relative Humidity:

Measurement Range: 0 - 100% Non-condensing

Accuracy: +/- 3 % to 100% RH

Resolution: 0.1% RH

### Data Retention:

Capacity: up to 48 hours per test

Sample rate: 1/min or 1/sec, automatic, parameter dependent

### Test Specifications:

#### Incubator Testing Specifications

- Measures temperature at 6 points as per IEC 60601-2-19 (sub-clause 50.102)
- Measures humidity at the center of the enclosure per IEC 60601-2-19 (sub-clause 103.1)
- Measures airflow at 4 points in the enclosure per IEC 60601-2-19 (sub-clause 104.1)
- Measures sound level inside the incubator per IEC 60601-2-19 (sub-clause 102.1)
- Measures sound level inside the incubator during alarm per IEC 60601-2-19 (sub-clause 102.2)
- Measures sound level outside the incubator per IEC 60601-2-19 (sub-clause 102.3)
- Measures incubator temperature at center point (A) of incubator as per IEC 60601-2-19 (sub-clause 50.102)
- Measures average incubator temperature (AIT) per IEC 60601-2-19 (sub-clause 2.9.106)
- Measures incubator temperature (IT) per IEC 60601-2-19 (sub-clause 2.9.105)
- Measures average temperature (AT) per IEC 60601-2-19 (sub-clause 2.9.103)
- Measures average temperature at 4 points per IEC 60601-2-19 (sub-clause 50.102)

- Measures skin temperature sensor accuracy per IEC 60601-2-19 (sub-clause 50.104)2
- Confirms skin temperature sensor control of incubator per IEC 60601-2-19 (sub-clause 50.105)
- Confirms independent incubator temperature sensor accuracy per IEC 60601-2-19 (sub-clause 50.106)
- Measures difference between average temperature and control temperature per IEC 60601-2-19 (sub-clause 50.107)
- Measures actual warm-up time compared to operating manual per IEC 60601-2-19 (sub-clause 50.108)
- Measures overshoot and recovery to steady state per IEC 60601-2-19 (sub-clause 50.109)
- Confirms humidity reading per IEC 60601-2-19 (sub-clause 50.110)

#### Radiant Warmer Testing Specifications<sup>1</sup>

- Measures skin temperature sensor accuracy per IEC 60601-2-21 (sub-clause 50.101)2
- Measures average warmer temperature at 4 points per IEC 60601-2-21 (sub-clause 50.102)
- Measures average warmer temperature at center point (A) of incubator as per IEC 60601-2-21 (sub-clause 50.102)
- Confirms skin temperature sensor control of warmer per IEC 60601-2-21 (sub-clause 50.103)
- Confirms airflow is less than 0.1 m/s per IEC 60601-2-21 (sub-clause 50.102)
- Measures sound level on the mattress per IEC 60601-2-21 (sub-clause 102.1)
- Measures sound level on the mattress during alarm per IEC 60601-2-21 (sub-clause 102.1)
- Measures sound level in front of the warmer per IEC 60601-2-21 (sub-clause 102.1)

<sup>1</sup> Radiant 1 warmer software to be released following general release of Incubtest

<sup>2</sup> Sensor accuracy heater is an optional accessory

*All specifications subject to change without notice.*



130 - 4020 Viking Way, Richmond BC, Canada V6V 2L4

Phone: 604.291.7747 • Fax: 604.294.2355

Toll-free (North America only): 800.667.6557

Email: customerservice@datrend.com

© DATREND SYSTEMS INC. MARCH 2015