



PSR-11-917-J

Response: 13 sec

Life: 60 mos

Warranty: 18 mos

OEM Equip

In vivo: 9445

Lifecare

PLV102

MSA:

Mini Ox I, II

Mini Ox 3, 3000

Puritan Bennett:

7820 Monitor

Respironics

Esprit

Precision Medical

PM5900 Monitor

SLE: T-801

Stephan: Christina

Taema:

Ceasar Vent

Eolo Vent

Teledyne:

AX 300 Analyzer

MX 300 Monitor



PSR-11-917-J1

Response: 13 sec

Life: 60 mos

Warranty: 18 mos

OEM Equip

Maxtec:

OM-25A Analyzer

OM-25E Monitor

OM-25MEL Monitor

MAX-250 ESF



PSR-11-917-M

Response: 13 sec

Life: 60 mos

Warranty: 18 mos

OEM Equip

Drager:

C2000 Incubator

Carefusion:

Avea Ventilator

SiPap Infant Vent

Infrasonics: Adult
Star Ventilator

Siemens:

300 / 900 Vents

Spacelabs Analyzer

Taema: Ceasar (old)

Versamed: iVent



PSR-11-917-M1

Response: 13 sec

Life: 60 mos

Warranty: 18 mos

OEM Equip

Datascope Monitor



PSR-11-917-M2

Response: 13 sec

Life: 60 mos

Warranty: 18 mos

OEM Equip

Puritan Bennett:

760 Ventilator

840 Ventilator



PSR-11-917-MH

Response: 6 sec

Life: 60 mos

Warranty: 18 mos

OEM Equip

Criticare: Poet

Hamilton Medical:

Arabella CPAP

Galileo Ventilator

Raphael Ventilator



PSR-11-917-MH1

Response: 6 sec

Life: 60 mos

Warranty: 14 mos

OEM Equip

Airshield: Incubator

ATOM Medical:

V-2200 Incubator

Drager:

C2000 Incubator

Advanced Oxygen Sensors for Respiratory

Proprietary Developments...

Enjoy competitive advantages
unmatched by competitors.

T90 Response...

Choose your response time without
giving up sensor life or increasing
dimensions.

Expected Life...

Our technology offers the longest
life and reliability in the industry.

Warranty...

A realistic warranty based on
the 60-70% average oxygen
level in normal service. Don't
be misled, get the facts.

In-House Production...

The widest range of sensors are
"Made in the USA".

Quality...

Our sensors are subject to the
most extensive stability testing
in the industry: output in air,
under 30" WC pressure and at
100% oxygen.

Service and Support...

99.5% of orders ship in 24 hours
and our phones are answered by
real live customer service reps
trained to help you.

**AII-11-60**

Response: 13 sec

Life: 60 mos

Warranty: 18 mos

OEM EquipAnalytical Industries

2000 Analyzer

2000A Analyzer

2000M Monitor

**AII-11-60HC**

Response: 13 sec

Life: 60 mos

Warranty: 18 mos

OEM EquipAnalytical Industries

2000HC Analyzer

**PSR-11-33**

Response: 30 sec

Life: 24 mos

Warranty: 12 mos

OEM EquipHamilton:

Aladdin

Amadeus

Veolar

**PSR-11-33-2**

Response: 30 sec

Life: 24 mos

Warranty: 12 mos

OEM EquipHudson:

5525 Analyzer

5550 Monitor

5552 Monitor

5828 Analyzer

**PSR-11-55**

Response: 9 sec

Life: 42 mos

Warranty: 14 mos

OEM EquipHudson:

5577 Analyzer

5590 Monitor

6477 Analyzer

**PSR-11-58**

Response: 13 sec

Life: 60 mos

Warranty: 18 mos

OEM EquipHudson:

5800 Monitor

5801 Analyzer

**PSR-11-58HC**

Response: 20 sec

Life: 60 mos

Warranty: 18 mos

OEM EquipHudson:

5802 Analyzer

**PSR-11-60-08**

Response: 6 sec

Life: 60 mos

Warranty: 18 mos

OEM EquipSensidyne:

1000E Analyzer

Datascope:

Minray AS 3000

**PSR-11-60-23**

Response: 6 sec

Life: 60 mos

Warranty: 18 mos

OEM EquipSensidyne:

5000M Monitor

**PSR-11-77**

Response: 6 sec

Life: 60 mos

Warranty: 18 mos

OEM EquipTeledyne:

60T Analyzer

191 Monitor

200 Monitor

**PSR-11-915-2i**

Response: 20 sec

Life: 60 mos

Warranty: 16 mos


OEM EquipDrager:

Baby Log Vent

I8000 Incubator

Evita 2 & 4 Vent

Evita XL

**PSR-11-915-G**

Response: 20 sec

Life: 60 mos

Warranty: 14 mos

OEM EquipOhmeda:

Giraffe Incubator

Quality ...

Our Quality Assurance program is independently certified and complies with the following regulatory standards:

U.S.: 510(k) No. K952736
ISO 9001:2008

Europe: Annex II of Medical Device Directive
93/42/EEC as amended by
2007/07/EC

Canada: ISO 13485:2003





PSR-11-75-KE1

Response: 13 sec

Life: 60 mos

Warranty: 18 mos

OEM Equip

Maxtec:

OM25A Analyzer

Vapor Therm:

Precision Flow



PSR-11-75-KE2

Response: 13 sec

Life: 60 mos

Warranty: 18 mos

OEM Equip

Puritan Bennett:

740 Vent (old)



PSR-11-75-KE4

Response: 13 sec

Life: 60 mos

Warranty: 18 mos

OEM Equip

Carefusion:

Vela Ventilator



PSR-11-75-KE6

Response: 13 sec

Life: 60 mos

Warranty: 18 mos

OEM Equip

SLE: 2000 Ventilator



PSR-11-75-KE7

Response: 13 sec

Life: 60 mos

Warranty: 18 mos

OEM Equip

Maxtec:

OM-25A Analyzer

OM25ME / MEL

MaxBlend

Newport Medical:

E360 Ventilator

E500 Ventilator



PSR-11-75-KE8

Response: 6 sec

Life: 60 mos

Warranty: 18 mos

OEM Equip

SLE: 5000 Ventilator



PSR-11-75-KE9

Response: 6 sec

Life: 60 mos

Warranty: 18 mos

OEM Equip

Carefusion:

Sentry Blender

Maxtec:

OM-25A Analyzer

OM-25ME / 25MEL

MaxBlend



PSR-11-75-KE10+

Response: 6 sec

Life: 60 mos

Warranty: 18 mos

OEM Equip

Maxtec

OMO2+A Analyzer

Distributed by:



Committed to Quality

100% QC Testing...

the most extensive testing in the industry provides you with the most reliable O2 sensors possible:

- Extended leak test confirms integrity of seals
- Mandated stabilization period to ensure consistency and reliability
- Pressure tested per ISO 7767
- Cycle from air to 100% to ensure stability and linearity
- Multiple uA and mV sensor output and stability tests during assembly

Advanced Technology...

- Longer life in traditional sizes
- Faster response time without giving up service life
- Proprietary processing designed to minimize internal contamination
- A design that simplifies assembly and enhances quality

Superior Service...

- Speak with real people, trained specifically to help you
- Knowledgeable engineers
- Delivery within 24 hours

Quality...

Our QA Program mandates that every warranty claim be thoroughly investigated and the customer furnished with a complete report.