

Electrocardiograph



JPD-E06 Six Channel ECG



Features

- Portable 6 channel ECG machine
- Working modes: Manual/Auto/Arrhythmia Analysis modes
- 112mm、3/6 channel format recording
- Built in rechargeable Li-ion battery
- Adapte 110-230V, 50/60Hz AC power supply
- Real-time ECG waveform, 16 patients save/replay
- Built in USB/RS232 interface
- Feasible paper size: 112mm roll/Z-fold
- 5.7" large LCD displays working parameters and waveform
- ECG archive software optional

Electrocardiograph



JPD-E06 Six Channel ECG

Specifications

Lead	Standard 12 leads	Sensitivity valve	≤20μV
Lead acquisition	12bit/1000Hz(synchronously 12 leads)	Sensitivity	5, 10, 20mm/mV
Work modes	MAN、AUTO、ANALYSIS Modes	Inter-channel interference	≤0.5mm
Filter	AC filter : 50Hz / 60Hz	Recording mode	3ch, 3ch+, 3ch+,, 6ch. Default Man. mode: 3ch. Default Auto. mode: 6ch
	EMG filter : 25Hz / 45Hz	Recording method	Thermal printing system
	Anti-drift filter : 0.15Hz(Adaptive)	Paper size	112mm roll/Z-fold
CMMR	≥120dB (with AC filter)	Paper speed	25、50mm/s(±3%)
Input circuit	Floating;	Case storage	16 ECG cases, expandable
	Protection circuit against defibrillator effect	Display	320x240 graphic 5.7 inch LCD
Input impedance	≥50MΩ	Dimension	330x270x75(mm)
Input circuit current	≤0.05μA	Weight	About 2.8Kg
Input dynamic range	±10mVp-p	Communication port	USB/RS232 port
Patient current leakage	<10μA	Fuse	T2A/250V, Φ5x20
Calibrating voltage	1mV±2%	Power supply	AC110~230V(±10%), 50/60Hz(±1Hz)
Voltage tolerance	±500mV		DC: built in rechargeable battery
Time constant	≥3.2s		14.4V(2000mAh)
Frequency response	0.05~160Hz(-3db)		
Noise level	<15μVp-p		

Configurations

Main unit	1 set
Patient cable	1 set
Limb electrodes	1 set(4 pcs)
Chest electrodes	1 set(6 pcs)
Power cable	1 set
Grounding cable	1 set
Adapter	1 pcs
Fuse	2 pcs
Paper axis	1 pcs
Printing paper	1 roll
Operation manual	1 copy



Shenzhen Jumper Medical Equipment Co., Ltd

※The material is subject to change without notice

5th Floor Building No. 34, Baoyuan Industrial Zone, Xixiang Street,
Baoan District, Shenzhen, Guangdong 518102 .P.R. China

Phone: +86-755-26696279 Fax: +86-755-26852025

Mail: info@jumper-medical.com