


# JPD-300P

FETAL MONITOR [ 7 inch ]



# JPD-300P

FETAL MONITOR [ 7 inch ]



## [ Product introduction ]

The monitor of lightweight and space-saving  
It can meet the need of whole stages of labor (including antepartum and intrapartum monitoring for fetus and mother)  
It also can be connected to the obstetrical central station

## [ Options ]

- ▲ Measurement of maternal SpO<sub>2</sub>,HR,NIBP
- ▲ Maternal temperature measurement
- ▲ Maternal ECG and PR measurement
- ▲ Twins monitor interface, Twins FHR measurement
- ▲ Build-in network/PC connecting interface
- ▲ FAS(fetal acoustic stimulator )

## [ Features ]

- ▲ Lightweight, space-saving, easy operation
- ▲ Interface of parameters of FHR , TOCO, FM, with digital and curve display
- ▲ Battery indication, pregnant No. display
- ▲ Alarm system
- ▲ 7-inch TFT color screen Folding up to 90 degree
- ▲ Build-in long life thermal printer
- ▲ DSP technology, realtime identify and measurement for FHR,accurate and reliable results
- ▲ Special high sensitive transducer
- ▲ 12 hours data storage which can be playback and printed
- ▲ Special patient event and clinical event marking bottom
- ▲ Low ultrasound power,safe to fetus
- ▲ Wide range of applying voltage(100-240V)
- ▲ Automatic Fetal Movement Detection



## SHENZHEN JUMPER MEDICAL EQUIPMENT CO.,LTD.

Add: 5th Floor,Building No.34, Baoyuan Industrial Zone,XixiangStreet, Baoan District,Shenzhen  
Tel: +86-755-26696279 Fax: +86-755-26852025 Email:info@jumper-medical.com  
Web: www.jumper-medical.com