

JPD-900W Blood Pressure Monitor

JUMPER

(Wrist type)



Features

JUMPER Wrist Full Automatic Blood Pressure Monitor is for use by medical professionals or at home and a non-invasive blood pressure measurement system intended to measure the diastolic, systolic blood pressures and pulse rate on adult individual by using a non-invasive technique in which an inflatable cuff is wrapped around the wrist. The cuff circumference is limited to 6.1023 inches to 9.8425 inches.

Packaging information

Unit dimension:

G/W(KG):0.19

N/W(KG):0.11

Size:108×90×90mm

Carton dimension: (48PCS/ carton)

G/W(KG):11

N/W(KG):9.5

Size:480×310×430mm